

JAY LowZone Gel Switch

Move with Confidence, Live with Freedom

The **JAY LowZone Gel Switch Cushion** combines a slow-moving gel core with temperature-sensitive visco-elastic memory foam for reliable pressure relief and all-day comfort. Its any-way-up Switch design allows flipping and rotating for even wear. Stable foam layers, orientation markings, and a protective inner liner enhance durability, while the water-resistant PU cover is easy to wipe clean or machine wash.



Product Type:
Gel and foam pressure-reducing wheelchair cushion

Construction:
Central slow-moving gel pouch encased in temperature-sensitive visco-elastic memory foam

Cover Type: Wipedown cover with internal Plaitition liner

Depth / Height 3" (8 cm)

Colour Black/White

Maximum User Weight 158kg

Size 17" x 17", 18" x 16", 20" x 18", 22" x 18"



Sales, Hire, Service & Repairs

t 08 7119 7777

www.thrivemobility.au

Unit 6, 220 North East Road, Klemzig 5087
klemzig@thrivemobility.au

See our YouTube Channel for product reviews & more